

GaitKeeperGo Quick Start Guide

Recording your first 4-metre walk test — follow the numbers around the photo

STAND FACING SUBJECT

Open GaitKeeperGo and tap "Start 4M Walk."



CHECK POSITION

Subject on/behind the green box. Pitch & Roll near 0°.



SCAN THE FLOOR

Scan the floor at the subject's feet to load the walkway.



PLACE THE WALKWAY

Tap to place it. A green box marks the start line.



TURN AND WALK 4M

Walk to the far end and stand on the red circle.



START RECORDING

Tap **Start**. Call "3-2-1-Walk" — subject walks at a normal pace. Follow on screen instructions.



STOP THE TEST

Tap **STOP** once the subject fully crosses the red line.



REVIEW THE RESULT

Tap "View Reports" for falls risk, speed trend & MobilityDNA.